

Volumetric Analysis/Titration Validation Test 2020

Time allowed:

minutes

45

Instructions

Please ensure you enter your name and circle your teacher's initials below.
Scientific calculators only. Chemistry

(Nam				
T	eacher:	(circle)			
\	CEM	NM	KL	MX	

Question 1 17 marks

A bottle of store bought cloudy ammonia was analysed to determine the concentration of ammonia present.

25.00 (\pm 0.03) mL of cloudy ammonia was pipetted into a 250 (\pm 0.15) mL volumetric flask and made up to the mark with distilled water.

A standard solution of hydrochloric acid was used, with a concentration of 0.103 (\pm 0.002) mol L^{-1} .

The standardised HC ℓ was titrated against 20.00 (± 0.03) mL aliquots of diluted cloudy ammonia, each with 3-5 drops of bromophenol blue added.

The average titre volume obtained was 16.23 (± 0.10) mL.

a)	Write an ionic equation for the reaction between the household ammonia hydrochloric acid.	and the (1 mark)
b)	The cloudy ammonia advertises its ammonia concentration as 15-20 g/L. Calculate the concentration of ammonia in the cloudy ammonia solution a determine whether it fits within this range. Note: Error analysis not require	

c) Complete the table below, calculating the % error in each measurement.

Measurement	Value	% error
25mL pipette	25.00 (± 0.03)	0.12
20 mL pipette	20.00 (± 0.03)	0.15
burette	16.23 (± 0.10)	
Volumetric flask	250 (± 0.15)	
[HCl]	0.103 (± 0.002)	

	0.103 (± 0.002)	[HCl]	[H
(1 mark			
n part b) to the correct number	express your answer calculated the box provided below.		d)
(3 marks)			
	dy ammonia:	oncentration of NH₃ in cloud	Concei
	g/L	±	
	ourette with distilled water only, von of ammonia? Fully explain yo		e)

f)		Bromophenol blue is an indicator that changes colour in the range of pH 3-4.6. With the use of an appropriate equation, explain why this indicator was used.
		(3 marks)

Question 2 14 marks

Acid 'X' is a weak, diprotic, organic acid found naturally in many fruits. It is often used in a concentrated form as a flavouring to give food a sour taste.

A student wanted to set up a titration to determine the concentration of 'Acid X' present in various samples of grapefruit juice to compare their quality and flavour.

The student researched that grapefruit juice can have a percentage by mass of 'Acid X' in the range of 2.5 - 5.0 %. The molar mass for Acid 'X' is 134 gmol^{-1} .


The average mass of 20.00 mL samples of grapefruit juice was found to be 20.80 g.

The student had standardised a solution of sodium hydroxide and determined its concentration to be 0.105 molL⁻¹.

a)	Using the information above, perform appropriate calculations to determine if t grapefruit juice needs to be diluted prior to the titration. Explain your decision.	he
	graperrait juice riceus to be unuted prior to the titration. Explain your decision.	(6 marks)

b) In the laboratory, the student had access to a 20 mL and 25 mL pipette, as well as volumetric flasks of 125 mL, 250 mL and 500 mL. Using this equipment, give the

- c) Describe an assumption the student needed to make when planning their analysis of the grapefruit juice. (1 mark)
- d) The pH was monitored throughout the titration and the following curve was obtained.



i) Explain why there are two equivalence points.

(1 mark)

ii)	Using the formula H_2X to denote Acid 'X', write an ionic equation to describe the reaction that has gone to completion at the first equivalence point (12.5 mL of NaOH added).
iii)	On the titration curve, draw the curve you would expect to obtain if the concentration of the standard NaOH solution was doubled and a total of 50 mL of was added. (3 marks)

Question 3 13 marks

Antacids are medications that treat heartburn, indigestion and acid reflux. The active ingredients are carbonate containing compounds. A standardised solution of hydrochloric acid (1.13 molL⁻¹) was used to determine the percentage by mass of magnesium carbonate in a common antacid tablet using a back titration procedure. The assumption was made that all the carbonate present was in the form of magnesium carbonate.

Step 1: A 1.67 g sample of the antacid tablet was

Step 2: 50.0 mL of the standardized hydrochloric acid was added to the conical flask;

Step 3: At the completion of the reaction, the mixture in the conical flask was added to a 500 mL volumetric flask and made up to the mark with distilled water;

Step 4: A 25.00 mL aliquot of the diluted mixture was titrated with 0.0530 molL⁻¹ sodium hydroxide, using phenolphthalein as the indicator.

The scientists performed a rough titration and then recorded the results of a further 4 titrations in the table below:

	Titration 1	Titration 2	Titration 3	Titration 4
Final volume (mL)	33.80	37.90	31.70	28.80
Initial volume (mL)	5.25	9.50	3.20	0.30
Volume titrated (mL)	28.55	28.40	28.50	28.50

_		
a)	Calculate the average volume of sodium hydroxide used in the titration sequence.	
	Vav =	(1 mark)
b)	Write an equation for the reaction of magnesium carbonate with hydrochloric acid in	n Step 2.
		(1 mark)
c)	Write an ionic equation for the reaction taking place in Step 4.	
		(1 mark)

Calculate the percentage by mass of magnesium carbonate in the antacid tablet.	(8 marks)
When conducting the titration above, give the solutions that you would use to following pieces of equipment with: (2 marks)	o wash the

Equipment	Solution to be rinsed with
burette	
500 mL volumetric flask	