


Q: 20mL of the sample solution is pipetted into a 250mL conical flask and diluted to 150mL with distilled water. It's titrated with 0.004mol L<sup>-1</sup> KIO<sub>3</sub> solution. The average titre is 7.65mL. Calculate the number of moles of ascorbic acid reacted. Calculate the concentration of ascorbic acid in mol L<sup>-1</sup>, mg/100mL and mg/100g.



Independent variable	Whether an open or unopened orange juice package is titrated.
Dependent variable	Concentration of vitamin C in the orange juice package.
Controlled variables	<ul style="list-style-type: none"> <li>• Surrounding temperature.</li> <li>• Surrounding pressure.</li> <li>• Volume of indicator added.</li> <li>• Indicator used.</li> </ul>

Systematic errors:

1. Inherent error in the burette.
2. Inherent error in the measuring cylinder.
3. Inherent error in the pipette.

(I)

Random errors:

1. Parallax error – viewing burette and measuring cylinders from different angles.
2. Endpoint determination – a visual endpoint is always slightly beyond the equivalence point due to needing to see the colour change.
3. Estimating values between graduations.

(PEE)

