Question Revision

1.

The concentration of hydrochloric acid in a liquid concrete cleaner was determined using the following method. A 25.00 mL sample was diluted to 250 mL in a volumetric flask. A pipette was used to transfer a 10.00 mL sample of this diluted concrete cleaner to a conical flask. An indicator was added and the sample of diluted concrete cleaner was titrated against a 0.200 M solution of sodium carbonate in a burette. The titre was found to be 23.15 mL.

2.

The contents of a laboratory reagent bottle containing solid potassium hydroxide, $KOH_{(s)}$, has been contaminated by moisture. You have been given the job of determining the extent of contamination. This is to be done by titrating a solution of this KOH with a standardised hydrochloric acid solution (HCl) of pH 1.14. You are to use the following procedure:

- Step 1: Dissolve 2.80 g of the contaminated KOH_(s) into 500.0 mL of solution.
- Step 2: Fill a burette with this solution and use it to titrate a 20.00 mL sample of the HCl_(aq) using a suitable indicator.
- **Step 3:** Repeat the titrations until a consistent end point is obtained.

After completing the titrations, the following data is available:

- Titration end point = 15.80 mL of $KOH_{(aq)}$.
- pH of $HCl_{(aq)} = 1.14$.

3.


A supermarket brand of vinegar is to be analysed for its ethanoic acid content by a titration technique. This requires the preparation of a Na₂CO₃(aq) primary standard. The primary standard is made by dissolving 1.416 g of anhydrous Na₂CO₃(s) in some distilled water and making the solution up to exactly 500.0 mL using a volumetric flask.

A 50.00 mL sample of the vinegar is diluted to exactly 1.000 L in another volumetric flask. Four 25.00 mL samples of this solution are placed into separate conical flasks and titrated with the $Na_2CO_3(aq)$ primary standard. On average 23.55 mL of the Na_2CO_3 solution was required for equivalence.

$$Na_2CO_3(aq) + 2CH_3COOH(aq) -----> 2NaCH_3COO(aq) + H_2O(l) + CO_2(g)$$

Determine the percentage by mass of ethanoic acid in the undiluted vinegar if this vinegar has a density of 1.060 g mL⁻¹.

- The concentration of hydrochloric acid in a liquid concrete cleaner was determined using the following method. A 25.00 mL sample was diluted to 250 mL in a volumetric flask. A pipette was used to transfer a 10.00 mL sample of this diluted concrete cleaner to a conical flask. An indicator was added and the sample of diluted concrete cleaner was titrated against a 0.200 M solution of sodium carbonate in a burette. The titre was found to be 23.15 mL.
- (a) Calculate the concentration of HCl in the original sample of concrete cleaner.



$$2HCl(aq) + Na_2CO_3(aq) \longrightarrow 2NaCl(aq) + H_2O(l) + CO_2(g)$$

$$\begin{split} &n(Na_2CO_3)_A=cV=0.200 \ x \ 0.02315=0.00463 \\ &n(HCl)_B=2 \ x \ n(Na_2CO_3)_A=2 \ x \ 0.00463=0.00926 \ mol \\ &c(HCl)_C=(n(HCl)_B)/0.0100=0.00926/0.01=0.926 \ M \\ &n(HCl)_C=cV=0.926 \ x \ 0.250=0.2315 \ mol \\ &c(HCl)_D=n/V=0.2315/0.025=\textbf{9.26} \ \textbf{M} \end{split}$$

<u>2.</u>


The contents of a laboratory reagent bottle containing solid potassium hydroxide, $KOH_{(s)}$, has been contaminated by moisture. You have been given the job of determining the extent of contamination. This is to be done by titrating a solution of this KOH with a standardised hydrochloric acid solution (HCl) of pH 1.14. You are to use the following procedure:

- Step 1: Dissolve 2.80 g of the contaminated KOH_(s) into 500.0 mL of solution.
- Step 2: Fill a burette with this solution and use it to titrate a 20.00 mL sample of the HCl_(aq) using a suitable indicator.
- Step 3: Repeat the titrations until a consistent end point is obtained.

After completing the titrations, the following data is available:

- Titration end point = 15.80 mL of $KOH_{(aq)}$.
- pH of $HCl_{(aq)} = 1.14$.

Calculate the percentage by mass of moisture in the contaminated KOH(s)



$$\begin{split} &n(HCl)_A = cV = 10^{\text{-}1.14} \text{ x } 0.020 = 1.45 \text{ x } 10^{\text{-}3} \text{ mol} \\ &n(KOH)_B = n(HCl)_A = 1.45 \text{ x } 10^{\text{-}3} \text{ mol} \\ &c(KOH)_B = n/V = (1.45 \text{ x } 10^{\text{-}3})/0.01580 = 9.17 \text{ x } 10^{\text{-}2} \text{ M} \\ &c(KOH)_C = 9.17 \text{ x } 10^{\text{-}2} \text{ M} \\ &n(KOH)_C = c \text{ x } V = (9.17 \text{ x } 10^{\text{-}2}) \text{ x } 0.500 = 0.0459 \\ &m(NaOH)_D = n \text{ x } M = 0.0459 \text{ x } 56.1056 = 2.58 \text{ g} \\ &\%(KOH) \text{ in impure sample} = (2.58/2.80) \text{ x } 100 = \textbf{92.0}\% \end{split}$$

3.


A supermarket brand of vinegar is to be analysed for its ethanoic acid content by a titration technique. This requires the preparation of a $Na_2CO_3(aq)$ primary standard. The primary standard is made by dissolving 1.416 g of anhydrous $Na_2CO_3(s)$ in some distilled water and making the solution up to exactly 500.0 mL using a volumetric flask.

A 50.00 mL sample of the vinegar is diluted to exactly 1.000 L in another volumetric flask. Four 25.00 mL samples of this solution are placed into separate conical flasks and titrated with the $Na_2CO_3(aq)$ primary standard. On average 23.55 mL of the Na_2CO_3 solution was required for equivalence.

$$Na_2CO_3(aq) + 2CH_3COOH(aq) -----> 2NaCH_3COO(aq) + H_2O(l) + CO_2(g)$$

Determine the percentage by mass of ethanoic acid in the undiluted vinegar if this vinegar has a density of 1.060 g mL⁻¹.

[10 marks]



```
\begin{split} &n(Na_2CO_3)_A = m/M = 1.416/105.9884 = 1.336 \text{ x } 10^{-2} \text{ mol} \\ &c(Na_2CO_3)_B = n/V = (1.336 \text{ x } 10^{-2})/0.500 = 2.67 \text{ x } 10^{-2} \text{ M} \\ &n(Na_2CO_3)_{added \text{ from Burette } C} = cV = (2.67 \text{ x } 10^{-2}) \text{ x } 0.02355 = 6.29 \text{ x } 10^{-4} \text{ mol} \\ &n(CH_3COOH)_D = 2 \text{ x } n(Na_2CO_3)_{added \text{ from Burette } C} = 2 \text{ x } (6.29 \text{ x } 10^{-4}) = 1.26 \text{ x } 10^{-3} \text{ mol} \\ &c(CH_3COOH)_E = n/V = (1.26 \text{ x } 10^{-3})/0.02500 = 5.03 \text{ x } 10^{-2} \text{ M} \\ &n(CH_3COOH)_E = cV = (5.03 \text{ x } 10^{-2}) \text{ x } 1 = 5.03 \text{ x } 10^{-2} \text{ mol} \\ &c(CH_3COOH)_E = n/V = (5.03 \text{ x } 10^{-2})/0.0500 = 1.01 \text{ M} \end{split} density of vinegar = 1060 g / L therefore, mass of 1 L solution = 1060 g
```

therefore, mass of 1 L solution = 1060 g number of moles CH_3COOH in 1 L = 1.01 therefore, mass of CH_3COOH in 1 L = n x M = 1.01 x 60.0519 = 60.5g %(CH_3COOH) = [(mass of CH_3COOH)/(mass of solution)] x 100 = (60.5/1060) x 100 = **5.7**%